
# FASTest<sup>®</sup> LH ad us. vet.

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	before LH peak (prooestrus)		LH peak (oestrus)	after LH peak (metoestrus)	
<b>LH test</b>	variably negative – positive if positive, test repetition after 2 h → negative		positive test repetition after 24 h → positive	negative test repetition after 2 h → negative	
<b>LH value ng/ml</b>	pulsatile LH emission fluctuates between < 1 and > 1		<b>praeovulatory LH &gt; 1</b>	LH basal values < 1	
<b>cycle state</b>	early prooestrus	late prooestrus	oestrus	early metoestrus	late metoestrus
<b>progesterone value ng/ml</b>	basal value: < 0,5	↑ up to ca. 1	↑ increase > 2-4	ovulation: 4-10	after ovulation: > 10
<b>vaginal cytology</b>	parabasal/basal/intermediate cells, esp. erythrocytes	mainly superficial cells and erythrocytes	> 90% superficial cells, „acaryote clods“, background “clear”	superficial cells ↓ ↑ parabasal/intermediate cells, partly neutrophil granulocytes	↑ parabasal/intermediate cells, superficial cells ↓ < 20%
<b>vaginoscopy</b>	rose mucosa / ↑ oedema sanguinary secretion	pinkish mucosa / ↑ oedema sanguinary secretion	pale mucosa/oedema maximal clear and less secretion	pinkish mucosa, hyperplastic, sticky yellowish mucous secretion	rose mucosa, moist, lucent marginal-no secretion
<b>bearing</b>	partial toleration of male dogs	toleration of male dogs	toleration of male dogs + mating	rejection of male dogs	