MegaELISA[®] CRP canine ad us. vet.

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In vitro Diagnostikum

GEBRAUCHSINFORMATION

Enzyme immunoassay for the quantitative detection
of C-reactive protein in plasma or serum of the dog

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Contents	Page
1. Introduction	14
2. Intended use	14
3. Test principle	14
4. Test-kit components	15
5. Stability and Storage	16
6. Information on the sample material	16
7. Test preparation	17
8. Test procedure	18
9. Test evaluation	19
10. Safety measures and precautions	21
11. Liability	22
12. Short protocol MegaELISA® CRP canine	23
Abbreviations	2

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1. INTRODUCTION

As response to tissue injury and inflammation, the body produces, among others, so-called acute-phase proteins (APP) in the liver. The C-reactive protein (CRP) is the most sensitive and fastest acute-phase protein of the dog, with low serum concentrations during normal homeostasis and a quick answer during and after beginning of inflammation. The function of CRP is, among others, stimulation of the unspecific immune defence by activation of phagocytes and the complement cascade.

The measurement of CRP in dogs is useful for the diagnostics and follow-up of inflammatory diseases. In healthy dogs, there is a low CRP level below 10 mg/l. A drastic increase in the CRP level occurs on inflammatory processes within few hours (Ø 24–28 h by a factor of 10 to 100) due to infections, immune disorders, neoplasia and traumata. As soon as the inflammation subsides through a successful therapy, the CRP level decreases rapidly within Ø 24–28 h.

Therefore CRP can be regarded as a real-time marker. Indicating inflammations, CRP is more sensitive than the determination of rectal temperature or leucocyte differential diagnostics. Another advantage of CRP determination is its independence to endogenous or iatrogenic disturbing factors (e. g. glucocorticoids, stress).

The **MegaELISA® CRP** canine is based on highly specific antibodies against canine CRP. It can be carried out with dog plasma or serum in the veterinary practice and/or in the laboratory.

2. INTENDED USE

MegaELISA® CRP canine is an enzyme immunoassay for the quantitative detection of C-reactive protein (CRP) in plasma or serum of the dog.

3. TEST PRINCIPLE

The **MegaELISA® CRP** canine is an indirect quantitative immunoassay to detect C-reactive protein (CRP) using a sandwich principle.

Specific polyclonal antibodies against CRP are bound to the surface of the wells. The diluted plasma or serum samples to be tested as well as the Controls are pipetted into the wells and incubated at room temperature (20–25°C). After a wash step, the Conjugate against CRP, labelled with horse radish peroxidase (HRP), is added and incubated at room temperature (20–25°C). Non-bound Conjugate is removed in another wash step. After addition of Substrate (Tetramethylbenzidine, TMB), with positive samples, the bound enzyme changes the colourless solution in the wells into a blue solution. By addition of Stop solution (HCl), the colour changes from blue to yellow.

With help of the Standard curve, the results of the samples are quantifiable.

4. TEST-KIT COMPONENTS

4.1. REAGENTS / COMPONENTS PROVIDED WITH THE TEST-KIT

The reagents of one package are sufficient for 96 determinations.

	96's kit	Content	
1.	12 strips with 8 wells	ELISA plate , coated with specific antibodies against canine CRP	
2.	100 ml	Sample dilution buffer (SDB, contains 0.06 % Bronidox), ready to use, transparent cap	
3.	30 ml	Wash buffer (50× buffer), transparent cap	
4.	12 ml	Assay buffer , ready to use, transparent cap	
5.	0.4 ml	Standard stock solution (500 ng/100 µl canine CRP), red cap	
6.	0.05 ml	Control 1 ("low Control"), yellow cap	
7.	0.05 ml	Control 2 ("high Control"), yellow cap	
8.	12 ml	Conjugate (peroxidase-conjugated CRP-specific antibodies, ready to use, transparent cap	
9.	12 ml	Substrate (3,3',5,5'-Tetramethylbenzidine [TMB]/H ₂ O ₂), ready to use, brown cap. Do not expose to sun light.	
10.	12 ml	Stop solution (1 M hydrochloric acid [HCl]), ready to use, transparent cap	
11.	1	Microplate lid	
12.	1	Evaluation sheet	
13.	1	Certificate of Analysis	
14.	1	Instructions for use	

4.2. REQUIRED, NON-PROVIDED REAGENTS / DEVICES

- Distilled or deionised water
- Disposable sample tubes (sample dilution)
- Laboratory pipettes (range between 5 and 1000 µl)
- Multichannel pipette (100 µl)
- Vortex mixer
- Filter paper/laboratory cloths
- Measuring cylinder (1000 ml)
- Stop watch
- ELISA plate shaker (orbital shaker, 500 rpm)
- Washing station for ELISA plates (350 µl)
- Photometer for ELISA plates (450 nm, reference filter 590–650 nm)
- Waste bin with a 0.5 % hypochlorite solution

5. STABILITY AND STORAGE

Store at
2–8°C

Expiry date
– see label

For veterinary use
only

Lot number

In vitro diagnosticum

Do not use test-kit
components from
different kits, lot
numbers or beyond
stated expiry date.

Follow instructions for
use precisely

- Microbial contamination must be avoided. After the expiration date, the guarantee of quality is not provided.
- When opening the sealed aluminium pouch with the ELISA plate, care must be taken that the reclosable clip closure is not damaged or even separated.
- A direct influence of light onto the colourless Substrate must be absolutely avoided to prevent degradation or colour change to blue by autooxidation. If the Substrate colour has changed to blue, it cannot be used any more (functional incapacity).

6. INFORMATION ON THE SAMPLE MATERIAL

- Use only dog plasma or serum samples.
- Centrifuge whole blood samples within 4 hours.
- Keep in mind that the sample material, as well as all used test-kit components, should have reached room temperature at the time of application.
- Mix the sample material well before use.
- Non-cooled (**20–25°C**), plasma or serum should be tested within 24 hours. At **2–8°C**, the plasma or serum can be stored up to 5 days. Plasma or serum samples can be aliquoted and permanently stored at minimum –20°C. Avoid repeated freezing and thawing.
- Some lipaemic samples can cause a background colour. In such cases, it is recommended to repeat the test with better quality samples.
- Do not heat inactivate.
- **ATTENTION:** Partially filled and/or insufficient mixed EDTA, Citrate or Heparin tubes could create interference effects resulting in varying test results.

7. TEST PREPARATION

Read instructions for use carefully before starting the test. For the reliability of the results, it is important to follow the instructions for use exactly as described. Carry out the test in the given order and without delay. Use clean disposable pipette tips for each pipetting step.

7.1. GENERAL

- Remove the strips from the aluminium pouches after they have reached room temperature. To avoid loss by evaporation, it is advisable to cover or mask them. Immediately return non-used strips into the pouch, close the pouch and store it at 2–8°C.
- Mix the reagents well immediately before their use.
- After use, reagents should be stored at 2–8°C again.
- Strips that were once used cannot be used again.
- Damaged test-kit components must not be used.
- Direct contact of samples with the test-kit components has to be avoided (cross-contamination).
- Direct solar radiation during testing must be avoided.

7.2. PREPARATION OF WASH BUFFER 1:50

Mix 1 part of the Wash buffer concentrate with 49 parts of distilled water. When crystals are visible in the Wash buffer concentrate, dissolve these by warming the Wash buffer concentrate in a water bath at 37°C. The diluted buffer will be stable for 4 weeks if stored at 2–8°C.

7.3. SAMPLE PREPARATION / SAMPLE DILUTION 1:100

- The samples to be tested are diluted 1:100 with Sample dilution buffer (SDB, e.g. 10 µl sample + 990 µl SDB) before beginning. The controls must also be diluted!
- Homogenise the sample-buffer mixture (SBM) well (repeated mixture with pipette or vortex mixer).
- **ATTENTION:** During the use of an ELISA washing station, the SBM must be free of particles. To guarantee this, please centrifuge the SBM at 5000 rpm for 5 minutes.
- To determine very low CRP concentrations (< 5 mg/l) in the sample and Control, the samples should be diluted 1:50 or 1:25. Lower dilutions than 1:25 are not recommended.
- To determine very high CRP concentrations (> 50 mg/l), the samples should be measured at higher dilutions (e.g., 1:200 or 1:500).

7.4. PREPARATION OF THE STANDARDS

- The standards have to be prepared freshly before each use. Use the SDB for this.
- The undiluted standard has a concentration of 500 ng/100 µl.
- Prepare several standard dilutions.

ID	Concentration	Standard dilution
Standard A	500 ng/ml	900 µl SDB + 100 µl Standard stock solution
Standard B	278 ng/ml	400 µl SDB + 500 µl Standard A
Standard C	154 ng/ml	400 µl SDB + 500 µl Standard B
Standard D	86 ng/ml	400 µl SDB + 500 µl Standard C
Standard E	48 ng/ml	400 µl SDB + 500 µl Standard D
Standard F	0 ng/ml	SDB

Table 1: Preparation of the standard dilutions

8. TEST PROCEDURE

- Mix sample material well before application.
- **It is recommended to test controls, standards and samples in duplicates.**
- Before beginning, the position of the patient samples, the Controls and the standard dilutions on the ELISA plate must be determined on the Evaluation sheet supplied with the kit.
- Use a multichannel pipette for Assay buffer, Conjugate, TMB and stop solution!

8.1. SAMPLING AND FIRST INCUBATION

- Place the required amount of strips into the holding frame.
- Add **100 µl** each of **Assay buffer** into each well (multichannel pipette).
- Add **20 µl** each of
 - **Standard A to F**
 - **diluted Controls**
 - **diluted samples (SBM)**
 in the desired order into the particular well.
- Cover the ELISA plate and incubate for **60 minutes at room temperature (20–25°C)** on the shaker (500 rpm).

8.2. WASHING – CAREFULLY and STRICTLY FOLLOWING THE INSTRUCTIONS!

AUTOMATIC WASHING, 3× with 350 µl each

- Remove the incubation solution by emptying into a waste bin with hypochlorite. Disposal must be made according to official regulations.
- **ATTENTION:** SBM still containing any particles should be manually removed from the cavities by shaking off (danger of blockage of the washing needles).

- Observe correct setting of the washing station (**3×350 µl**).
- To reach optimal washing results use at least 350 µl Wash buffer per well and wash step.
- Aspirate liquid completely after each wash step.
- After the final wash step, beat the plate carefully but thoroughly on filter paper/laboratory cloth. There should not be any residual moisture in the cavities.

8.3. SECOND INCUBATION

- Add **100 µl Conjugate** into each well (multichannel pipette).
- Cover the ELISA plate and incubate for **30 minutes at room temperature (20–25 °C)** and **in the dark** (e. g., by covering with aluminium foil) on the shaker (500 rpm).

8.4. WASHING

- Wash as described in issue 8.2., but with **5 washing steps**.

8.5. THIRD INCUBATION

- Add **100 µl Substrate (TMB)** into each well (multichannel pipette).
- Cover the ELISA plate and incubate for **15–20 minutes at room temperature (20–25 °C)** and **in the dark** on the shaker (500 rpm).
- Add **100 µl Stop solution** into each well, in the same order and at the same rate as for the TMB substrate solution (multichannel pipette).

8.6. MEASUREMENT

- After careful mixing (tipping softly at the plate rim), the extinction can be measured at 450 nm within 10 minutes after addition of the Stop solution. An additional measurement at the reference wavelength 590–650 nm is recommended.
- If duplicates or multiple determinations were done, calculate the **mean extinction values**.

9. TEST EVALUATION

9.1. STANDARD CURVE

The standard curve is created by plotting the standard concentration on the x-axis (linear scale) against the extinction of the standards on the y-axis (linear scale) (see fig.1). A 4-parameter method (curve fitting) should be used for automatic data reduction. Alternatively, a quadratic (polynomial) analysis is also possible.

ID	Concentration	Extinction at 450 nm
Standard A	500 ng/ml	3.088
Standard B	278 ng/ml	1.297
Standard C	154 ng/ml	0.456
Standard D	86 ng/ml	0.211
Standard E	48 ng/ml	0.114
Standard F	0 ng/ml	0.006

Table 2: Example for evaluation of the standard concentration


Fig.1: Example for evaluation of the standard concentration

ID	Value
A	0.032593
B	1.85576
C	849.366
D	11.2598
R ²	0.9996

Calculation by means of curve fitting:

$$y = (A - D) / (1 + (x/C)^B) + D$$

$$y = (0.032593 - 11.2598) / (1 + x/849.366)^{1.85576} + 11.2598$$

9.2. INTERNAL QUALITY CONTROL

The values of the controls are shown on the CoA supplied. These should lie on the determined standard curve.

A variation of the required values as well as any clouding or blue staining of the colourless Substrate before addition into the cavities can be a hint on expiry of reagents. If the prescribed issues are not fulfilled, following should be ensured before test repetition:

- Shelf life of the used reagents
- Functionality of the used machines (calibration)
- Visual control of the test-kit components for contamination or leakage; a blueish stained Substrate solution should not be used.

9.3. CUT-OFF

- CRP in healthy dogs < 10 mg/l
- Range of testing 5–50 mg/l
- Sensitivity: 2.9 mg/l

9.4. EVALUATION OF TEST RESULTS

The CRP concentration of the dog samples is calculated by multiplication of the value read from the standard curve by the dilution used for the respective sample. Values > 10 mg/l are considered elevated/acute inflammation.

Samples with a higher absorbance than standard A should be retested with a higher dilution.

9.5. LIMITATIONS OF THE TEST PROCEDURE

- The diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptoms as well as serological data.
- Bacterial contamination or repeated freeze-thaw cycles of the specimen may affect the extinction values.

9.6. INTERFERENCES

There were no cross reactions with mouse, goat, horse, rat, cattle, chicken, sheep, pig and human.

There were no interferences with hemoglobin/hemolysate.

10. SAFETY MEASURES AND PRECAUTIONS

- Wear safety clothes (disposable gloves, lab coat, safety glasses) and wash the hands after finishing the test procedure.
- No eating, drinking and smoking during test procedure.
- The sample material and Controls must be seen as potentially infectious and disposed of accordingly after test procedure, together with the used test-kit components.
- Do not use reagents of other manufacturers together with the reagents of this test-kit.
- Do not mix reagents of different Lot numbers.
- Do not change the caps of the individual reagents among each other.
- Close bottles immediately after use to avoid evaporation and microbial contamination.
- After first opening, the reagent bottles must be checked for microbial contamination before reusing.

- – Sample dilution buffer contains 0.06 % Bronidox as preservative. Avoid contact with skin, eyes and clothing. Upon contact, wash with water. If ingested, call a physician.
- – Stop solution contains 1 M Hydrochloric acid. Hydrochloric acid is caustic and can cause severe burns. Upon contact with eyes or skin, rinse thoroughly with water and consult a doctor. Keep out of the reach of children.
- – Substrate solution contains hydrogen peroxide.

Absolutely avoid eye, skin or mucosa contact.

- **ATTENTION:** Liquid waste containing Stop solution must be neutralised prior to addition to a hypochlorite solution.
- Only use clean pipette tips, dispenser and laboratory material.
- Never pipet with your mouth.
- Label sample material and associated sample tube to ensure a precise assignment.
- Use a new sample tube for each sample.
- During the incubation and wash steps, the cavities must be completely covered with liquid.
- The cavities of the ELISA plate must not stay empty longer as necessarily required.
- To avoid cross contamination and falsely elevated test results, pipet the patient samples and conjugate carefully into the cavities of the ELISA plate.
- **MegaELISA® CRP** canine is only provided for the use of qualified staff that knows the work technique objectionable.

The test procedure, the information, the safety measures and warnings in the Instructions for use are to be followed strictly. For use of the test-kit in diagnostic instruments, the test method has to be validated. Each change in appearance, composition and the test procedure as well as each use in combination with other products not authorised by the manufacturer is not permitted; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents being associated herewith. No liability is accepted for false test results due to visual evaluation.

11. LIABILITY

The entire risk of liability in connection with the use of this product remains with the purchaser. The manufacturer is not liable for indirect, special or consequential damages of any kind resulting from the usage, test performance and evaluation of this product.

12. SHORT PROTOCOL MegaELISA® CRP canine

TEST PREPARATION

Warm needed strips and reagents to **room temperature (20–25°C)**.

Dilute Wash buffer 1:50 with distilled water.

Dilute samples and controls 1:100 with Sample dilution buffer
(e. g. 10 μl sample + 990 μl SDB)

Prepare standard dilution series

TEST PROCEDURE

Add **100 μl Assay buffer**
into each well.

Add **20 μl of 1:100 diluted plasma or serum samples** or **Controls**
or **Standard dilutions**
into each well.

Cover and incubate **60 minutes at room temperature (20–25°C) on the shaker.**

Wash 3 \times with 350 μl Wash buffer each.

Add **100 μl Conjugate**
into each well.

Cover and incubate **30 minutes at room temperature (20–25°C)**
in the dark on the shaker.

Wash 5 \times with 350 μl 1 \times Wash buffer each.

Add **100 μl Substrate (TMB)**
into each well.

Cover and incubate **15–20 minutes at room temperature (20–25°C)**
in the dark on the shaker.

Add **100 μl Stop solution**
into each well.

Photometric evaluation of the cavities at 450 nm.

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