# MegaELISA<sup>®</sup> EHRLICHIA canis ad us. vet.

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*In vitro* Diagnostikum

## GEBRAUCHSINFORMATION

Enzyme immunoassay for the qualitative detection  
of IgG antibodies against *Ehrlichia canis*  
in plasma or serum of the dog

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Contents	Page
1. Introduction	13
2. Intended use	13
3. Test principle	13
4. Test-kit components	14
5. Stability and Storage	15
6. Information on the sample material	15
7. Test preparation	16
8. Test procedure	17
9. Test evaluation	18
10. Features of performance	19
11. Safety measures and precautions	20
12. Liability	21
13. Short protocol <b>MegaELISA® EHRlichia canis</b>	22
Abbreviations	23

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## 1. INTRODUCTION

Canine monocytic ehrlichiosis (CME) is caused by the rickettsia *Ehrlichia canis*, which are mainly transmitted by the brown dog tick (*Rhipicephalus sanguineus*). *Ehrlichia canis* is found in many parts of the world, especially in the Mediterranean area, but also in Switzerland and partly Germany. Ehrlichiosis is common in the dog, but seldom in humans.

The concentration of specific antibodies increases sharply 14 to 21 days post infection. A four-fold titre increase of antibodies in an interval of two weeks (seroconversion) is indicative for an acute infection. In the acute stage, the dog shows apathy, anorexia, fever and lymphadenitis. In the subclinical phase, clinical symptoms are missing, but not the typical ehrlichiosis laboratory results like hyperglobulinaemia and thrombocytopenia. Chronic phase animals show slight up to life-threatening symptoms: spontaneous bleedings, neurological disorders, anaemia, severe loss of weight as well as spleno- and hepatomegaly. Indirect antibody detection is known to be an important diagnostic tool diagnosing CME beside clinical symptoms, case history (travel abroad) and direct antigen detection.

**MegaELISA® EHRlichIA canis** is based on highly specific antigens for the fast and reliable detection of antibodies against *Ehrlichia canis* in plasma or serum of infected dogs.

---

## 2. INTENDED USE

**MegaELISA® EHRlichIA canis** is an enzyme immunoassay for the qualitative detection of IgG antibodies against *Ehrlichia canis* in plasma or serum of the dog.

---

## 3. TEST PRINCIPLE

The **MegaELISA® EHRlichIA canis** is a qualitative immunoassay to detect IgG antibodies against *Ehrlichia canis* using a sandwich principle.

Specific *Ehrlichia canis* antigens are bound to the surface of the well in the ELISA plate. A dilution of the dog plasma or serum to be tested as well as controls are incubated in a first step at room temperature (20–25°C) in the wells of the ELISA plate. After a wash step, anti-dog IgG conjugate labelled with horse radish peroxidase (HRP) is added and incubated at room temperature (20–25°C). Non-bound conjugate is removed in another wash step. After addition of Substrate (Tetramethylbenzidine [TMB]), with positive samples, the bound enzyme changes the colourless solution in the wells of the ELISA plate into a blue solution. By addition of Stop solution (Sulfuric acid), the colour changes from blue to yellow. The extinction is proportional to the concentration of *Ehrlichia*-specific IgG antibodies present in the sample.

## 4. TEST-KIT COMPONENTS

### 4.1. COMPONENTS / REAGENTS PROVIDED WITH THE TEST-KIT

The reagents of one package are sufficient for 48 or 96 determinations.

Content	48's kit	96's kit
1. <b>ELISA plate</b> with separable test strips in holding frame with 8 determinations each, coated with <i>Ehrlichia canis</i> antigens	6 separable strips	12 separable strips
2. <b>Sample dilution buffer</b> (SDB, pH 7.2 ± 0.2, contains 0.045 % Proclin® 300), ready to use, <b>blue cap</b>	25 ml	50 ml
3. <b>Wash buffer</b> (10× TRIS buffer, pH 7.2 ± 0.2, contains Tween® 20, 0.045 % Proclin® 300), <b>white cap</b>	25 ml	50 ml
4. <b>Positive Control</b> (contains 0.045 % Proclin® 300), ready to use, <b>red cap</b>	0.5 ml	1 ml
5. <b>Negative Control</b> (contains 0.045 % Proclin® 300), ready to use, <b>blue cap</b>	0.5 ml	1 ml
6. <b>Cut-off Control</b> (contains 0.045 % Proclin® 300), ready to use, <b>green cap</b>	0.5 ml	1 ml
7. <b>Conjugate</b> (peroxidase-conjugated anti-dog IgG antibodies, pH 7.2 ± 0.2, contains 0.045 % Proclin® 300), ready to use, <b>white cap</b>	6 ml	12 ml
8. <b>Substrate</b> (3,3',5,5'-tetramethylbenzidine [TMB]/ H <sub>2</sub> O <sub>2</sub> ), ready to use, <b>yellow cap. Do not expose to sun light!</b>	6 ml	12 ml
9. <b>Stop solution</b> (0.2 M Sulfuric acid), ready to use, <b>red cap</b>	8 ml	15 ml
10. <b>Cover foil</b>	1	1
11. <b>Evaluation sheet</b>	1	1
12. <b>Instructions for use</b>	1	1

### 4.2. REQUIRED, NON-PROVIDED REAGENTS / DEVICES

- Distilled or deionised water
- Disposable sample tubes
- Laboratory pipettes (5–100 µl, 1 ml)
- Vortex mixer
- Filter paper/laboratory cloths
- Measuring cylinder (1000 ml)
- Stop watch

- Washing station for ELISA plates or multichannel pipette (300 µl)
- Photometer for ELISA plates (450 nm, reference filter 620 nm)
- Waste bin with a 0.5 % hypochlorite solution

## 5. STABILITY AND STORAGE

Store at  
2–8 °C

Expiry date  
– see label

For veterinary use  
only

Lot number

For research only

Do not use test-kit  
components from  
different kits, lot  
numbers or beyond  
stated expiry date.

Follow instructions for  
use precisely

Microbial contamination must be avoided. After the expiration date, the guarantee of quality is not provided.

When opening the sealed aluminium pouch with the ELISA plate, care must be taken that the reclosable clip closure is not damaged or even separated! Remove the required amount of ELISA strips, leave the remaining strips in the bag, close the bag and store it at 2–8 °C.

A direct influence of light onto the colourless Substrate must be absolutely avoided to prevent degradation/autooxidation/colour change to blue. If the Substrate colour has changed to blue, it cannot be used any more (functional incapacity).

## 6. INFORMATION ON THE SAMPLE MATERIAL

- Use only dog plasma or serum samples!
- Mix the sample material well before use!
- Non-cooled (**20–25 °C**), plasma or serum should be tested within 4 hours! At **2–8 °C**, the plasma or serum can be stored up to 5 days. Plasma or serum samples can be aliquoted and permanently stored at minimum –20 °C. Avoid repeated freezing and thawing.
- Do not heat inactivate!
- Keep in mind that the sample material, as well as all used test-kit components, should have reached room temperature at the time of application.
- **ATTENTION:** Partially filled and/or insufficient mixed EDTA, Citrate or Heparin tubes could create interference effects resulting in varying test results.

## 7. TEST PREPARATION

**!** Read instructions for use carefully before starting the test. For the reliability of the results, it is important to follow the instructions for use exactly as described. Carry out the test in the given order and without delay. Use clean disposable tips for each pipetting step.

---

### 7.1. GENERAL

- Remove the test strips from the pouches after they have reached room temperature. To avoid loss by evaporation, it is advisable to cover or mask them.
- Mix the reagents well immediately before their use.
- After use, both non-used test strips (in closed aluminium pouch) and reagents should be stored at 2–8°C again.
- Test strips that were once used cannot be used again.
- Damaged test-kit components must not be used.
- Direct contact of samples with the test-kit components has to be avoided (cross-contamination).
- Direct solar radiation during testing must be avoided.

### 7.2. PREPARATION OF WASH BUFFER 1:10

Mix 1 part of the Wash buffer concentrate with 9 parts of distilled water. When crystals are visible in the Wash buffer concentrate, dissolve these by warming the Wash buffer concentrate up to 37°C in a water bath. The diluted buffer will be stable for 5 days if stored at room temperature (20–25°C).

### 7.3. SAMPLE PREPARATION / SAMPLE DILUTION 1:101 OR 1+100

- Plasma and serum has to be diluted for testing 1:101 or 1+100 with SDB, e. g. 5 µl sample + 500 µl SDB.
- Homogenise the sample-buffer mixture (SBM) well (repeated mixture with pipette or vortex mixer).
- **ATTENTION:** During the use of an ELISA washing station, the SBM must be free of particles! To guarantee this, please centrifuge the SBM with 5000 rpm for 5 minutes.
- Controls are ready for use and need not to be diluted.

## 8. TEST PROCEDURE

- Mix sample material well before application!
- Before beginning, the position of the patient samples and the standards/controls in the microtitre plate must be determined on the result sheet supplied with the kit. **It is recommended to test Controls and samples in duplicates.**
- Use clean disposable pipette tips for every pipetting step.

### 8.1. SAMPLING AND FIRST INCUBATION

- Place the required amount of test strips into the holding frame.
- Add **50 µl** each of
  - **Positive Control**
  - **Negative Control**
  - **Cut-off Control (2×)**
  - **diluted samples (SBM)**in the desired order into the particular cavity of the test strip. **Leave one well empty for the Substrate Blank.**
- Cover the ELISA plate and incubate for **60 minutes at room temperature (20–25°C)**.

### 8.2. WASHING – CAREFULLY and STRICTLY FOLLOWING THE INSTRUCTIONS!

#### MANUAL WASHING

- Remove the incubation solution by emptying into a waste bin with hypochlorite. Disposal must be made according to official regulations!
- Beat the ELISA plate carefully several times onto filter paper/laboratory cloth.
- Wash 4× with 300 µl Wash buffer per cavity. Shake plate with Wash buffer lightly, then beat it.
- **Drain completely between each wash step, then beat carefully several times on a dry unused spot!**
- Before the last beating take care that the reagent for the next pipetting step is prepared.

#### AUTOMATIC WASHING

- **ATTENTION:** SBM still containing any particles should be manually removed from the cavities by shaking off (danger of blockage of the washing needles).
- Observe correct setting of the washing station!
- To reach optimal washing results use at least 300 µl Wash buffer per cavity and wash step.
- Aspire liquid completely after each wash step.
- After the final wash step, beat the plate carefully but thoroughly on filter paper/laboratory cloth! There should not be any residual moisture in the cavities!

A variation of the required values as well as any clouding or blue staining of the colourless Substrate before addition into the cavities can be a hint on expiry of reagents. If the prescribed values are not fulfilled, following should be ensured before test repetition:

- Stability of the used reagents
- Functionality of the used machines (calibration)
- Visual control of the test-kit components for contamination or leakage; a blueish stained Substrate solution should not be used.

## 9.2. CALCULATION OF TEST RESULTS

The Cut-off value is calculated as the mean value of the two determined Cut-off controls.

**Cut-off = (Extinction Cut-off Control 1 + Extinction Cut-off Control 2) / 2**

e. g.  $(0.23 + 0.19) / 2 = 0.42 / 2 = 0.21$ ; Cut-off = 0.21

## 9.3. EVALUATION OF TEST RESULTS

### Positive test result

Sample with an extinction more than 10 % above the calculated Cut-off.

### Marginal test result

Sample with an extinction in the field of 10 % above and below the calculated Cut-off. A repetitive test after 2–4 weeks with a new sample being marginal, too, should be evaluated as negative!

### Negative test result

Sample with an extinction more than 10 % below the calculated Cut-off.

*Results in MEGACOR Units (MU)*

$\frac{\text{Patient (mean) extinction value} \times 10}{\text{Cut-off}} = [\text{MEGACOR Units} = \text{MU}]$

Example:  $\frac{0.956 \times 10}{0.21} = 45.5 \text{ MU}$

Cut-off:	10 MU	Negative:	< 9 MU
Grey zone:	9–11 MU	Positive:	> 11 MU

## 10. FEATURES OF PERFORMANCE

### 10.1. LIMITATIONS OF THE TEST PROCEDURE

- Bacterial contamination or repeated freeze-thaw cycles of the specimen may affect the extinction values.
- Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data.

## 11. SAFETY MEASURES AND PRECAUTIONS

- Wear safety clothes (disposable gloves) and wash the hands after finishing the test procedure.
  - No eating and drinking during test procedure.
  - The sample material and Controls must be seen as potentially infectious and disposed of accordingly, together with the used test-kit components, after the test procedure.
  - Do not use reagents of other manufacturers together with the reagents of this test-kit.
  - Do not mix reagents of different Lot numbers.
  - Do not change the caps of the individual reagents among each other.
  - Close bottles immediately after use to avoid evaporation and microbial contamination.
  - After first opening, the reagent bottles must be checked for microbial contamination before reusing.
    - Sample dilution buffer, Wash buffer, Controls and Substrate contain < 0,0014 % mixture of CMIT/MIT as preservative.
    - Stop solution contains 0.2 M Sulfuric acid. Sulfuric acid irritates eyes and skin. Upon contact with the eyes, rinse thoroughly with water and consult a doctor! Keep out of the reach of children.
    - Substrate solution contains hydrogen peroxide.
- Absolutely avoid eye, skin or mucosa contact!**
- **ATTENTION:** Liquid waste containing Stop solution must be neutralised prior to addition to a hypochlorite solution!
  - Only use clean pipette tips, dispenser and laboratory material.
  - Never pipet with your mouth!
  - Label sample material and associated sample tube to ensure a precise assignment.
  - Use a new sample tube for each sample.
  - During the incubation and wash steps, the cavities must be completely covered with liquid.
  - The cavities of the ELISA plate must not stay empty longer as necessarily required.
  - To avoid cross contamination and falsely elevated test results, pipet the patient samples and conjugate carefully into the cavities of the ELISA plate.
  - **MegaELISA® EHRlichia** canis is only provided for the use of qualified staff that knows the work technique objectionable.

The test procedure, the information, the safety measures and warnings in the Instructions for use are to be followed strictly. For use of the test-kit in diagnostical instruments, the test method has to be validated. Each change in appearance, composition and the application as well as each use in combination with other products not authorised by the manufacturer is not permitted. The manufacturer is not liable for false results and incidents being associated herewith. The user himself is responsible for such changes. No liability is accepted for false rest results due to visual evaluation.

---

## **12. LIABILITY**

**The entire risk due to the performance of this product is assumed by the purchaser. The manufacturer shall not be liable for indirect, special or consequential damages of any kind resulting from the use of this product.**

### 13. SHORT PROTOCOL MegaELISA® EHRlichia canis

#### TEST PREPARATION

Warm needed test strips and reagents to **room temperature (20–25°C)**.

**Dilute Wash buffer 1:10** with distilled water.

**Dilute samples 1:101** with Sample dilution buffer (e.g. 5 µl sample + 500 µl SDB)

---

#### TEST PROCEDURE

Add **50 µl of Positive, Negative and 2× Cut-off Controls**  
or **1:101 diluted plasma or serum samples**  
each in a cavity.

Cover and incubate **60 minutes at room temperature (20–25°C)**.

**Wash 4×** with soft shaking with **300 µl Wash buffer each**.

Add **100 µl of Conjugate**  
into each cavity. Mix carefully.

Cover and incubate **45 minutes at room temperature (20–25°C)**.

**Wash 4×** with soft shaking with **300 µl 1× Wash buffer each**.

Add **100 µl of Substrate (TMB)**  
in each cavity. Mix carefully.

Cover and incubate **exactly 10 minutes at room temperature (20–25°C)**  
**in the dark**.

Add **100 µl of Stop solution**  
in each cavity.

**Photometric evaluation of the cavities at 450 nm.**

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