# MegaELISA® ENCEPHALITOZOON cuniculi ad us. vet.

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*In vitro* Diagnostikum

## GEBRAUCHSINFORMATION

Enzyme immunoassay for the qualitative detection  
of IgG antibodies against *Encephalitozoon cuniculi*  
in plasma or serum of the rabbit

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## Contents

	Page
1. Introduction	13
2. Intended use	13
3. Test principle	13
4. Test-kit components	14
5. Stability and Storage	15
6. Information on the sample material	15
7. Test preparation	16
8. Test procedure	17
9. Test evaluation	18
10. Safety measures and precautions	19
11. Liability	20
12. Short protocol <b>MegaELISA® ENCEPHALITOOZON cuniculi</b>	21
Abbreviations	2

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## 1. INTRODUCTION

The infection with *Encephalitozoon cuniculi* occurs via spores in the urine of infected animals. The uptake of the extremely resistant spores mostly occurs perorally by contaminated food or inhalation in contaminated surroundings. That is why the spreading in house and laboratory animals is significantly faster than in wild animals. Carnivores also could be infected by infected prey.

Neurological symptoms, signs of kidney failure, and eye skin inflammation may occur due to lesions in the central nervous system, kidney, and eye. Diseased rabbits can experience one or more of these symptoms. The most common is the so-called vestibular syndrome (head tilt, disorders of movement coordination and eye tremors). Other neurological symptoms can include seizures, incomplete paralysis, and loss of balance. In rare cases, increased aggression and loss of hearing or vision may also occur. In some animals, renal insufficiency with rather unspecific symptoms (loss of appetite or weight, dehydration, disorders of the mineral balance and bone metabolism as well as apathy) can be indicative.

The *in vivo* diagnosis of encephalitozoonosis in rabbits is problematic because of the large number of animals with chronic, asymptomatic infection. In this case, antibodies can be detected in the blood for years. The serological detection of antibodies against *E. cuniculi* is usually the most sensitive method to determine an infection; a negative antibody titre during this period normally excludes an infection with *E. cuniculi*.

**MegaELISA® ENCEPHALITOOZON cuniculi** is based on highly specific antigens for the fast and reliable detection of IgG antibodies against *Encephalitozoon cuniculi* in plasma or serum of infected rabbits.

---

## 2. INTENDED USE

The **MegaELISA® ENCEPHALITOOZON cuniculi** is an enzyme immunoassay for the qualitative detection of IgG antibodies against *Encephalitozoon cuniculi* in plasma or serum of the rabbit.

---

## 3. TEST PRINCIPLE

The **MegaELISA® ENCEPHALITOOZON cuniculi** is an indirect qualitative immunoassay to detect IgG antibodies against *Encephalitozoon cuniculi*.

Specific *Encephalitozoon cuniculi* antigens are bound to the surface of the wells. A dilution of the rabbit plasma or serum to be tested as well as controls are incubated in a first step at room temperature (20–25 °C) in the wells. After a wash step, anti-rabbit IgG conjugate labelled with horse radish peroxidase (HRP) is added and incubated at room temperature (20–25 °C). Non-bound conjugate is removed in another wash step. After addition of Substrate (Tetramethylbenzidine, TMB), with positive samples, the bound enzyme changes the colourless solution in the wells into a blue solution. By addition of Stop solution (Sulfuric acid), the colour changes from blue to yellow. The extinction is proportional to the concentration of *Encephalitozoon*-specific IgG antibodies present in the sample.

## 4. TEST-KIT COMPONENTS

### 4.1. REAGENTS / COMPONENTS PROVIDED WITH THE TEST-KIT

The reagents of one package are sufficient for 48 or 96 determinations.

	48's kit	96's kit	Content
1.	6 strips with 8 wells	12 strips with 8 wells	<b>ELISA plate</b> , coated with <i>Encephalitozoon cuniculi</i> antigens
2.	25 ml	50 ml	<b>Sample dilution buffer</b> (SDB, pH 7.2 ± 0.2, contains 0.045 % ProClin® 300), ready to use, <b>blue cap</b>
3.	25 ml	50 ml	<b>Wash buffer</b> (10× TRIS buffer, pH 7.2 ± 0.2, contains 0.045 % ProClin® 300), <b>white cap</b>
4.	0.5 ml	1 ml	<b>Positive Control</b> (contains 0.045 % ProClin® 300), ready to use, <b>red cap</b>
5.	0.5 ml	1 ml	<b>Negative Control</b> (contains 0.045 % ProClin® 300), ready to use, <b>blue cap</b>
6.	6 ml	12 ml	<b>Conjugate</b> (peroxidase-conjugated anti-rabbit IgG antibodies, pH 7.2 ± 0.2, contains 0.045 % ProClin® 300), ready to use, <b>white cap</b>
7.	6 ml	12 ml	<b>Substrate</b> (3,3',5,5'-tetramethylbenzidine [TMB]/ H <sub>2</sub> O <sub>2</sub> ), ready to use, <b>yellow cap. Do not expose to sun light.</b>
8.	8 ml	15 ml	<b>Stop solution</b> (0.2 M Sulfuric acid), ready to use, <b>red cap</b>
9.	1	1	<b>Cover foil</b>
10.	1	1	<b>Evaluation sheet</b>
11.	1	1	<b>Instructions for use</b>

### 4.2. REQUIRED, NON-PROVIDED REAGENTS / DEVICES

- Distilled or deionised water
- Disposable sample tubes
- Laboratory pipettes (range between 5 and 1000 µl)
- Vortex mixer
- Filter paper/laboratory cloths
- Measuring cylinder (1000 ml)
- Stop watch
- Washing station for ELISA plates or multichannel pipette (300 µl)
- Photometer for ELISA plates (450 nm, reference filter 620 nm)
- Waste bin with a 0.5 % hypochlorite solution

## 5. STABILITY AND STORAGE

Store at  
2–8°C

Expiry date  
– see label

For veterinary use  
only

Lot number

*In vitro* diagnosticum

Do not use test-kit  
components from  
different kits, lot  
numbers or beyond  
stated expiry date.

Follow instructions for  
use precisely

- Microbial contamination must be avoided. After the expiration date, the guarantee of quality is not provided.
- When opening the sealed aluminium pouch with the ELISA plate, care must be taken that the reclosable clip closure is not damaged or even separated.
- A direct influence of light onto the colourless Substrate must be absolutely avoided to prevent degradation or colour change to blue by autooxidation. If the Substrate colour has changed to blue, it cannot be used any more (functional incapacity).

## 6. INFORMATION ON THE SAMPLE MATERIAL

- Use only rabbit plasma or serum samples.
- Keep in mind that the sample material, as well as all used test-kit components, should have reached room temperature at the time of application.
- Mix the sample material well before use.
- Original samples, that is, undiluted plasma or serum, should be tested non-cooled (20–25°C) within 4 hours. At 2–8°C, the undiluted plasma or serum can be stored up to 5 days. Plasma or serum samples can be aliquoted and permanently stored at minimum –20°C. Avoid repeated freezing and thawing.
- Plasma or serum, diluted 1+600 in sample dilution buffer SDB (according to “7.3. sample preparation”, p. 16), can be stored for testings in **MegaELISA® ENCEPHALITOOZON cuniculi**. The sample-buffer mixtures (SBM) are stable for up to 28 days at room temperature, at 2–8°C they are stable up to 2 months, at –20°C they can be stored permanently and thawed and frozen up to 10 times.

**ATTENTION:** This information refers to the use of samples exclusively in the **MegaELISA® ENCEPHALITOOZON cuniculi** and the stability of the *Encephalitozoon* IgG antibodies in the SBM, this does not apply to other parameters.

- Some lipaemic or haemolytic samples can cause a background colour. In such cases, it is recommended to repeat the test with better quality samples.

- Do not heat inactivate.
- **ATTENTION:** Partially filled and/or insufficient mixed EDTA, Citrate or Heparin tubes could create interference effects resulting in varying test results.

---

## 7. TEST PREPARATION

**Read instructions for use carefully before starting the test. For the reliability of the results, it is important to follow the instructions for use exactly as described.**

**Carry out the test in the given order and without delay. Use clean disposable tips for each pipetting step.**

---

### 7.1. GENERAL

- Remove the required wells from the aluminium pouches after they have reached room temperature. To avoid loss by evaporation, it is advisable to cover or mask them. Immediately return non-used wells into the pouch, close the pouch and store it at 2–8°C.
- Mix the reagents well immediately before their use.
- After use, reagents should be stored at 2–8°C again.
- Wells that were once used cannot be used again.
- Damaged test-kit components must not be used.
- Direct contact of samples with the test-kit components has to be avoided (cross-contamination).
- Direct solar radiation during testing must be avoided.

### 7.2. PREPARATION OF WASH BUFFER 1:10

Mix 1 part of the Wash buffer concentrate with 9 parts of distilled water. When crystals are visible in the Wash buffer concentrate, dissolve these by warming the Wash buffer concentrate in a water bath at 37°C. The diluted buffer will be stable for 5 days if stored at room temperature (20–25°C).

### 7.3. SAMPLE PREPARATION / SAMPLE DILUTION 1:600

- Plasma and serum has to be diluted for testing 1:600 with SDB, please use at least 5 µl sample. We recommend to work with a pre-dilution.

Example:

Pre-dilution (1:10): 5 µl sample + 45 µl SDB.

Final dilution (1:600): 5 µl pre-dilution + 295 µl SDB (= final sample-buffer mixture SBM)

- Homogenise the SBM well (repeated mixture with pipette or vortex mixer).
- **ATTENTION:** During the use of an ELISA washing station, the SBM must be free of particles. To guarantee this, please centrifuge the SBM at 5000 rpm for 5 minutes.
- Controls are ready to use and need not to be diluted.

## 8. TEST PROCEDURE

- Mix sample material well before application.
- **It is recommended to test controls and samples in duplicates.**
- Before beginning, the position of the patient samples and the Controls in the wells must be determined on the Evaluation sheet supplied with the kit.

### 8.1. SAMPLING AND FIRST INCUBATION

- Place the required amount of wells into the holding frame.
- Add **50 µl** each of
  - **Positive Control**
  - **Negative Control**
  - **diluted samples (SBM)**in the desired order into the particular well. **Leave one well empty for the Substrate Blank.**
- Cover the ELISA plate and incubate for **60 minutes at room temperature (20–25 °C).**

### 8.2. WASHING – CAREFULLY AND STRICTLY FOLLOWING THE INSTRUCTIONS!

#### MANUAL WASHING

- Remove the incubation solution by emptying into a waste bin with hypochlorite. Disposal must be made according to official regulations.
- Beat the ELISA plate carefully several times onto filter paper/laboratory cloth.
- Wash 5× with 300 µl Wash buffer per well. Shake plate with Wash buffer slightly, then beat it.
- **Drain completely between each wash step, then beat carefully several times on a dry unused spot.**
- Before the last beating, take care that the reagent for the next pipetting step is prepared.

#### AUTOMATIC WASHING

- **ATTENTION:** SBM still containing any particles should be manually removed from the wells by shaking off (danger of blockage of the washing needles).
- Observe correct setting of the washing station (5×300 µl).
- To reach optimal washing results use at least 300 µl Wash buffer per well and wash step.
- Aspire liquid completely after each wash step.
- After the final wash step, beat the plate carefully but thoroughly on filter paper/laboratory cloth. There should not be any residual moisture in the wells.

### 8.3. SECOND INCUBATION

- Add **100 µl Conjugate** into each well except for the Substrate Blank well.
- Cover the ELISA plate and incubate for **45 minutes at room temperature (20–25 °C).**

### 8.4. WASHING

- Wash as described in issue 8.2.

### 8.5. THIRD INCUBATION

- Add **100 µl Substrate (TMB)** into each well.
- Cover the ELISA plate and incubate for exactly 15 minutes at room temperature (20–25°C) and **in the dark**.
- Add **100 µl Stop solution** into each well, in the same order and at the same rate as for the TMB substrate solution.

### 8.6. MEASUREMENT

- After careful mixing (tipping softly at the plate rim), the extinction can be measured at 450 nm within 30 minutes after addition of the Stop solution. An additional measurement at the reference wavelength 620 nm is recommended.
- Adjust the ELISA Plate Reader to zero using the Substrate Blank. If, due to technical reasons, this is not possible, subtract the extinction value of this well from all other extinction values measured in order to obtain reliable results.
- If duplicates or multiple determinations were done, calculate the **mean extinction values**.
- **ATTENTION:** Strongly positive samples can cause blackish precipitates of the Substrate solution. Pre-dilution of the sample with physiological NaCl solution, for example 1+1, is recommended. Then dilute the sample 1:600 with Sample dilution buffer and multiply the extinction values by 2.

---

## 9. TEST EVALUATION

### 9.1. REAGENT STABILITY – INTERNAL QUALITY CONTROL

During each testing, Positive and Negative Controls should be performed to ensure the stability of reagents as well as the correct test procedure.

The test was correct, when following reference values were reached:

- **Substrate Blank**      **Extinction < 0.1**
- **Negative Control**    **Extinction < 0.1**
- **Positive Control**     **0.8 < Extinction < 2.2**

A variation of the required values as well as any clouding or blue staining of the colourless Substrate before addition into the wells can be a hint on expiry of reagents. If the prescribed issues are not fulfilled, following should be ensured before test repetition:

- Shelf life of the used reagents
- Functionality of the used machines (calibration)
- Visual control of the test-kit components for contamination or leakage; a blueish stained Substrate solution should not be used.

### 9.2. CUT-OFF

The Cut-off value is between 0.200 and 0.500.

### 9.3. EVALUATION OF TEST RESULTS

#### Positive test result

Sample with an extinction above the Cut-off of 0.5.

#### Marginal test result

Sample with an extinction in the field between 0.2 and 0.5 (Cut-off). A repetitive test after 2–4 weeks with a new sample being marginal, too, should be evaluated as negative.

#### Negative test result

Sample with an extinction below the Cut-off of 0.2.

### 9.4. INTERPRETATION OF TEST RESULTS

- The diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptoms as well as serological data.
- Bacterial contamination or repeated freeze-thaw cycles of the specimen may affect the extinction values.

---

## 10. SAFETY MEASURES AND PRECAUTIONS

- Wear safety clothes (disposable gloves, lab coat, safety glasses) and wash the hands after finishing the test procedure.
  - No eating, drinking and smoking during test procedure.
  - The sample material and Controls must be seen as potentially infectious and disposed of accordingly after the test procedure, together with the used test-kit components, after the test procedure.
  - Do not use reagents of other manufacturers together with the reagents of this test-kit.
  - Do not mix reagents of different Lot numbers.
  - Do not change the caps of the individual reagents among each other.
  - Close bottles immediately after use to avoid evaporation and microbial contamination.
  - After first opening, the reagent bottles must be checked for microbial contamination before reusing.
    - Sample dilution buffer, Wash buffer, Controls and Substrate contain 0.045 % Pro-Clin® 300 as preservative.
    - Stop solution contains 0.2 M Sulfuric acid. Sulfuric acid irritates eyes and skin. Upon contact with the eyes, rinse thoroughly with water and consult a doctor. Keep out of the reach of children.
    - Substrate solution contains hydrogen peroxide.
- Absolutely avoid eye, skin or mucosa contact.**
- **ATTENTION:** Liquid waste containing Stop solution must be neutralised prior to addition to a hypochlorite solution.

- Only use clean pipette tips, dispenser and laboratory material.
- Never pipet with your mouth.
- Label sample material and associated sample tube to ensure a precise assignment.
- Use a new sample tube for each sample.
- During the incubation and wash steps, the wells must be completely covered with liquid.
- The wells must not stay empty longer as necessarily required.
- To avoid cross contamination and falsely elevated test results, pipet the patient samples and conjugate carefully into the wells.
- **MegaELISA® ENCEPHALITOOZON** cuniculi is only provided for the use of qualified staff that knows the work technique objectionable.

The test procedure, the information, the safety measures and warnings in the Instructions for use are to be followed strictly. For use of the test-kit in diagnostical instruments, the test method has to be validated. Each change in appearance, composition and the test procedure as well as each use in combination with other products not authorised by the manufacturer is not permitted; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents being associated herewith. No liability is accepted for false test results due to visual evaluation.

---

## 11. LIABILITY

**The entire risk of liability in connection with the use of this product remains with the purchaser. The manufacturer is not liable for indirect, special or consequential damages of any kind resulting from the usage, test performance and evaluation of this product.**

## 12. SHORT PROTOCOL MegaELISA® ENCEPHALITOOZON cuniculi

### TEST PREPARATION

Warm needed wells and reagents to **room temperature (20–25°C)**.

**Dilute Wash buffer 1:10** with distilled water.

**Dilute samples 1:600** with Sample dilution buffer.

---

### TEST PROCEDURE

Add **50 µl of Positive and Negative Controls**  
or **1:600 diluted plasma or serum samples**  
each into a well.

▼  
Cover and incubate **60 minutes at room temperature (20–25°C)**.

▼  
**Wash 5×** with soft shaking with **300 µl Wash buffer each**.

▼  
Add **100 µl of Conjugate**  
into each well except Substrate Blank. Mix carefully.

▼  
Cover and incubate **45 minutes at room temperature (20–25°C)**.

▼  
**Wash 5×** with soft shaking with **300 µl 1× Wash buffer each**.

▼  
Add **100 µl of Substrate (TMB)**  
into each well, also the **Substrate Blank**. Mix carefully.

▼  
Cover and incubate **exactly 15 minutes at room temperature (20–25°C)**  
**in the dark**.

▼  
Add **100 µl of Stop solution**  
into each well.

▼  
**Photometric evaluation of the wells at 450 nm.**  
Reference wavelength 620 nm.

▼  
Adjustment against Substrate Blank if possible.

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