MegaELISA[®] FeLV Antigen ad us. vet.

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In vitro Diagnostikum

GEBRAUCHSINFORMATION

Enzyme immunoassay for the qualitative detection
of Feline Leukaemia Virus (FeLV) antigens in plasma
or serum of the cat

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Contents	Page
1. Introduction	13
2. Intended use	13
3. Test principle	13
4. Test-kit components	14
5. Stability and storage	15
6. Information on the sample material	15
7. Test preparation	16
8. Test procedure	16
9. Test evaluation	18
10. Safety measures and precautions	20
11. Liability	21
12. Short protocol MegaELISA® FeLV Antigen	22
Abbreviations	2

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1. INTRODUCTION

Feline Leukaemia Virus (FeLV) is a highly contagious oncogenic RNA virus that causes both neoplastic and non-neoplastic diseases in cats world-wide.

The infection with Feline Leukaemia Viruses is most of all transmitted by secretions (e.g. saliva), excretions (e.g. feces, urine) as well as intrauterine by blood transfusion and via colostrum. As Feline Leukaemia Viruses can survive in the environment only for minutes, successful contamination only can happen during close contact to infected cats.

The course of an FeLV infection essentially depends on immune status, age as well as infection dose and virulence of the virus. Diseases caused by FeLV include lymphosarcoma, myelogenous leukaemia, thymic atrophy, non-regenerative anaemia and panleukopenia-like disease.

Because FeLV is immunosuppressive, it predisposes infected cats to a variety of secondary diseases. FeLV infected, viraemic cats normally show high concentrations of extracellular (free) p27 antigen 3 weeks post infection. p27 antigen detection is proved as the method of choice in FeLV diagnostics.

Based on highly specific monoclonal antibodies against p27 antigen, **MegaELISA® FeLV Antigen** is an important diagnostic tool for the evaluation of clinical suspicious cats (separation of FeLV carrier cats) and for FeLV check prior to vaccination.

Neither maternal antibodies nor a FeLV vaccination influence the test result of the MegaELISA® FeLV Antigen.

2. INTENDED USE

MegaELISA® FeLV Antigen is an enzyme immunoassay for the qualitative detection of Feline Leukaemia Virus (FeLV) in plasma or serum of the cat.

3. TEST PRINCIPLE

The **MegaELISA® FeLV Antigen** is an indirect qualitative immunoassay to detect Feline Leukaemia Virus (FeLV) antigen using a sandwich principle.

Monoclonal antibodies against p27 antigen are bound to the surface of the wells. The plasma or serum samples to be tested as well as the Controls and the Conjugate against p27, labelled with horse radish peroxidase (HRP), are pipetted into the wells and incubated at room temperature (20–25°C). Non-bound serum and Conjugate are removed by a wash step. After addition of Substrate (Tetramethylbenzidine, TMB), with positive samples, the bound enzyme changes the colourless solution in the wells into a blue solution. By addition of Stop solution (Sulfuric acid), the colour changes from blue to yellow. The extinction is proportional to the concentration of the FeLV antigens present in the sample.

4. TEST-KIT COMPONENTS

4.1. REAGENTS / COMPONENTS PROVIDED WITH THE TEST-KIT

The reagents of one package are sufficient for 48 or 96 determinations.

	48's kit	96's kit	Content
1.	6 strips with 8 wells	12 strips with 8 wells	ELISA plate , coated with monoclonal antibodies against p27. The wells of the individual strips can be broken off individually as required.
2.	12.5 ml	25 ml	Wash buffer (10× TRIS buffer, pH 7.2 ± 0.2, contains 0.045 % ProClin® 300), white cap
3.	0.25 ml	0.5 ml	Positive Control (contains 0.045 % ProClin® 300), ready to use, red cap
4.	0.25 ml	0.5 ml	Negative Control (contains 0.045 % ProClin® 300), ready to use, blue cap
5.	3.5 ml	7 ml	Conjugate (peroxidase-conjugated anti-p27 antibodies, contains 0.045 % ProClin® 300), ready to use, white cap
6.	6 ml	12 ml	Substrate (3,3',5,5'-tetramethylbenzidine [TMB]/H ₂ O ₂), ready to use, yellow cap. Do not expose to sun light.
7.	8 ml	15 ml	Stop solution (0.2 M Sulfuric acid), ready to use, red cap
8.	1	1	Cover foil
9.	1	1	Evaluation sheet
10.	1	1	Instructions for use

4.2. REQUIRED, NON-PROVIDED REAGENTS / DEVICES

- Distilled or deionised water
- Disposable sample tubes
- Laboratory pipettes (range between 5 and 1000 µl)
- Vortex mixer
- Filter paper/laboratory cloths
- Measuring cylinder (1000 ml)
- Stop watch
- Washing station for ELISA plates or multichannel pipette (300 µl)
- Photometer for ELISA plates (450 nm, reference filter 620 nm)
- Waste bin with a 0.5 % hypochlorite solution

5. STABILITY AND STORAGE

Store at
2–8°C

Expiry date
– see label

For veterinary use
only

Lot number

In vitro diagnosticum

Do not use test-kit components
from different kits, lot numbers
or beyond stated expiry date.

Follow instructions for
use precisely

- Microbial contamination must be avoided. After the expiration date, the guarantee of quality is not provided.
- When opening the sealed aluminium pouch with the ELISA plate, care must be taken that the reclosable clip closure is not damaged or even separated.
- A direct influence of light onto the colourless Substrate must be absolutely avoided to prevent degradation or colour change to blue by autooxidation. If the Substrate colour has changed to blue, it cannot be used any more (functional incapacity).

6. INFORMATION ON THE SAMPLE MATERIAL

- Use only cat plasma or serum samples.
- Keep in mind that the sample material, as well as all used test-kit components, should have reached room temperature at the time of application.
- Mix the sample material well before use.
- Non-cooled (**20–25°C**), plasma or serum should be tested within 4 hours. At **2–8°C**, the plasma or serum can be stored up to 5 days. Plasma or serum samples can be aliquoted and permanently stored at minimum **–20°C**. Avoid repeated freezing and thawing.
- Some lipaemic or haemolytic samples can cause a background colour. In such cases, it is recommended to repeat the test with better quality samples.
- Do not heat inactivate.
- **ATTENTION:** Partially filled and/or insufficient mixed EDTA, Citrate or Heparin tubes could create interference effects resulting in varying test results.

7. TEST PREPARATION

- !** Read instructions for use carefully before starting the test. For the reliability of the results, it is important to follow the instructions for use exactly as described. Carry out the test in the given order and without delay. Use clean disposable tips for each pipetting step.
-

7.1. GENERAL

- Remove the required wells from the aluminium pouches after they have reached room temperature. To avoid loss by evaporation, it is advisable to cover or mask them. Immediately return non-required wells into the pouch, close the pouch and store it at 2–8°C.
- Mix the reagents well immediately before their use.
- After use, reagents should be stored at 2–8°C again.
- Wells that were once used cannot be used again.
- Damaged test-kit components must not be used.
- Direct contact of samples with the test-kit components has to be avoided (cross-contamination).
- Direct solar radiation during testing must be avoided.

7.2. PREPARATION OF WASH BUFFER 1:10

Mix 1 part of the Wash buffer concentrate with 9 parts of distilled water. When crystals are visible in the Wash buffer concentrate, dissolve these by warming the Wash buffer concentrate in a water bath at 37°C. The diluted buffer will be stable for 5 days if stored at room temperature (15–25°C).

7.3. SAMPLE AND CONTROL PREPARATION

- Samples are applied undiluted.
 - **ATTENTION:** During the use of an ELISA washing station, the samples must be free of particles. To guarantee this, centrifuge the samples at 5000 rpm for 5 minutes.
 - Controls are ready to use and need not to be diluted.
-

8. TEST PROCEDURE

- Mix sample material well before application.
- **It is recommended to test controls and samples in duplicates.**
- Before beginning, the position of the patient samples and the Controls on the ELISA plate must be determined on the Evaluation sheet supplied with the kit.

8.1. SAMPLING AND FIRST INCUBATION

- Place the required amount of wells into the holding frame.
- Add **15 µl** each of
 - **Positive Control**
 - **Negative Control**
 - **undiluted sample**in the desired order into the particular well. **Leave one well empty for the Substrate Blank.**
- Add **50 µl Conjugate** into each well **except for the Substrate Blank.**
- Cover the ELISA plate and incubate for exactly **15 minutes at room temperature (20–25 °C).**

8.2. WASHING – CAREFULLY and STRICTLY FOLLOWING THE INSTRUCTIONS!

MANUAL WASHING

- Remove the incubation solution by emptying into a waste bin with hypochlorite. Disposal must be made according to official regulations.
- Beat the ELISA plate carefully several times onto filter paper/laboratory cloth.
- Wash 4× with 300 µl Wash buffer per well. Shake plate with wash buffer slightly, then beat it.
- **Drain completely between each wash step, then beat carefully several times on a dry unused spot.**
- Before the last beating, take care that the reagent for the next pipetting step is prepared.

AUTOMATIC WASHING

- **ATTENTION:** Samples still containing any particles should be manually removed from the wells by shaking off (danger of blockage of the washing needles).
- Observe correct setting of the washing station (4×300 µl).
- To reach optimal washing results use at least 300 µl Wash buffer per well and wash step.
- Aspire liquid completely after each wash step.
- After the final wash step, beat the plate carefully but thoroughly on filter paper/laboratory cloth. There should not be any residual moisture in the wells.

8.3. SECOND INCUBATION

- Add **100 µl Substrate (TMB)** into each well.
- Cover the ELISA plate and incubate for exactly **10 minutes at room temperature (20–25 °C) and in the dark.**
- Add **100 µl Stop solution** into each well, in the same order and at the same rate as for the TMB substrate solution.

8.4. MEASUREMENT

- After careful mixing (tipping softly at the plate rim), the extinction can be measured at 450 nm within 30 minutes after addition of the Stop solution. An additional measurement at the reference wavelength 620 nm is recommended.
- Adjust the ELISA Plate Reader to zero using the Substrate Blank. If, due to technical reasons, this is not possible, subtract the extinction value of this well from all other extinction values measured in order to obtain reliable results.
- If duplicates or multiple determinations were done, calculate the **mean extinction values**.
- **ATTENTION:** Strongly positive samples can cause blackish precipitates of the Substrate solution. Pre-dilution of the sample with physiological NaCl solution, for example 1+1, is recommended. Then multiply the extinction values by 2.

9. TEST EVALUATION

9.1. REAGENT STABILITY – INTERNAL QUALITY CONTROL

During each testing, Positive and Negative Controls should be performed to ensure the stability of reagents as well as the correct test procedure.

The test was correct, when following reference values were reached:

- **Substrate Blank** **Extinction < 0.1**
- **Negative Control** **Extinction < 0.1**
- **Positive Control** **0.7 < Extinction < 1.8**

A variation of the required values as well as any clouding or blue staining of the colourless Substrate before adding into the wells can be a hint on expiry of reagents. If the prescribed issues are not fulfilled, following should be ensured before test repetition:

- Shelf life of the used reagents
- Functionality of the used machines (calibration)
- Visual control of the test-kit components for contamination or leakage; a blueish stained Substrate solution should not be used.

9.2. CALCULATION OF THE TEST RESULTS

The test results of the samples (in percent) are calculated in relation to the Positive Control (PC), which corresponds to 100%.

$$\frac{\text{Extinction sample} \times 100}{\text{Extinction PC}} = xx \%$$

Example: Positive Control Extinction 1.2

Sample X Extinction 0.9

Sample Y Extinction 0.03

percentage calculation: Sample X: $\frac{0.9 \times 100}{1.2} = 75 \%$ Sample Y: $\frac{0.03 \times 100}{1.2} = 2.5 \%$

9.3. EVALUATION OF TEST RESULTS

Positive test result

Sample whose percentage value is above 10%.

Marginal test result

Sample whose percentage value is between 4 and 10%. A repetitive test after 2–4 weeks with a new sample being marginal, too, should be evaluated as negative.

Negative test result

Sample whose percentage value is below 4%.

9.4. INTERPRETATION OF THE TEST RESULTS

- The diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptoms as well as serological data.
- Bacterial contamination or repeated freeze-thaw cycles of the specimen may affect the extinction values.

MegaELISA® FeLV Antigen = NEGATIVE → no viraemia

- non-infected (approx. 95 %) or latent infected cat (no antigen detectable)
- test in the first 4 weeks post infection

MegaELISA® FeLV Antigen = 1st time POSITIVE → suspicion of viraemia

- transiently or persistently infected cat
- non-infected cat (false positive)

ADVICE: test repetition after 4–8 weeks. Isolation of the cat during this time is recommended (FeLV shedding!).

MegaELISA® FeLV Antigen = 2nd time POSITIVE → viraemia

- distinction transient or persistent viraemia
- 3rd test after another 6 weeks
- 4th test after another 10 weeks
- cat still positive → suspicion for persistent viraemia
 - **“progressor cat”** with high risk developing FeLV associated diseases
- cat becomes negative → suspicion of transient viraemia
 - **“regressor cat”** (complete virus elimination, healthy)
 - **latent infection** (integrated provirus in the bone marrow)

10. SAFETY MEASURES AND PRECAUTIONS

- Wear safety clothes (disposable gloves, lab coat, safety glasses) and wash the hands after finishing the test procedure.
 - No eating, drinking and smoking during test procedure.
 - The sample material and Controls must be seen as potentially infectious and disposed of accordingly after test procedure, together with the used test-kit components.
 - Do not use reagents of other manufacturers together with the reagents of this test-kit.
 - Do not mix reagents of different Lot numbers.
 - Do not change the caps of the individual reagents among each other.
 - Close bottles immediately after use to avoid evaporation and microbial contamination.
 - After first opening, the reagent bottles must be checked for microbial contamination before reusing.
 - – Wash buffer, Controls and Conjugate contain 0.045% Proclin® 300 as preservative.
 - – Stop solution contains 0.2 M Sulfuric acid. Sulfuric acid irritates eyes and skin. Upon contact with the eyes, rinse thoroughly with water and consult a doctor. Keep out of the reach of children.
 - – Substrate solution contains hydrogen peroxide.
- Absolutely avoid eye, skin or mucosa contact.**
- **ATTENTION:** Liquid waste containing Stop solution must be neutralised prior to addition to a hypochlorite solution.
 - Only use clean pipette tips, dispenser and laboratory material.
 - Never pipet with your mouth.
 - Label sample material and associated sample tube to ensure a precise assignment.
 - Use a new sample tube for each sample.
 - During the incubation and wash steps, the wells must be completely covered with liquid.
 - The wells must not stay empty longer as necessarily required.
 - To avoid cross contamination and falsely elevated test results, pipet the patient samples and conjugate carefully into the wells.
 - **MegaELISA® FeLV Antigen** is only provided for the use of qualified staff that knows the work technique objectionable.

The test procedure, the information, the safety measures and warnings in the Instructions for use are to be followed strictly. For use of the test-kit in diagnostical instruments, the test method has to be validated. Each change in appearance, composition and the test procedure as well as each use in combination with other products not authorised by the manufacturer is not permitted; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents being associated herewith. No liability is accepted for false test results due to visual evaluation.

11. LIABILITY

The entire risk of liability in connection with the use of this product remains with the purchaser. The manufacturer is not liable for indirect, special or consequential damages of any kind resulting from the usage, test performance and evaluation of this product.

12. SHORT PROTOCOL MegaELISA® FeLV Antigen

TEST PREPARATION

Warm needed wells and reagents to **room temperature (20–25°C)**.

Dilute Wash buffer 1:10 with distilled water.

The samples are used undiluted.

TEST PROCEDURE

Add **15 µl of Positive and Negative Controls**
or **15 µl of undiluted plasma or serum samples**
each into a well, except Substrate Blank.

↓
Add **50 µl Conjugate each**
into the wells, except Substrate Blank. Mix carefully.

↓
Cover and incubate **exactly 15 minutes at room temperature (20–25°C)**.

↓
Wash 4× with soft shaking with **300 µl 1× Wash buffer each**.

↓
Add **100 µl Substrate (TMB)**
in each well, **also the Substrate Blank**. Mix carefully.

↓
Cover and incubate **exactly 10 minutes at room temperature (20–25°C)**
in the dark.

↓
Add **100 µl Stop solution**
in each well.

↓
Photometric evaluation of the wells at 450 nm.
Reference wavelength 620 nm.

↓
Adjustment against Substrate Blank if possible.

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