MegaELISA® GIARDIA ad us. vet.

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In vitro Diagnostikum

GEBRAUCHSINFORMATION

Enzyme immunoassay for the qualitative detection
of *Giardia duodenalis* in feces of pocket pets, pets and farm animals

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Content	page
1. Introduction	13
2. Intended use	13
3. Test principle	14
4. Test-kit components	14
5. Stability and storage	15
6. Information on the sample material	16
7. Test preparation	16
8. Test procedure	17
9. Test evaluation	19
10. Safety measures and precautions	20
11. Liability	21
12. Short protocol MegaELISA® GIARDIA	22
Abbreviations	2

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Optional: Referenzwellenlänge 620 nm.

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1. INTRODUCTION

Giardia is known to be one of the most common enteritic parasites in pocket pets, pets, farm and wild animals as well as in humans (zoonosis) world-wide. The species *Giardia duodenalis* occurs in the different animal species and humans in varying genotypes (genotypes A to G), which differ in their infection and their host spectrum. Types A and B have zoonotic potential, whereas the others are more or less host-specific. Newborns and young animals are mostly affected. Prevalences vary in cats and dogs, depending on age (> 70% under 1 year), husbandry (10% in single husbandry up to 100% in breedings and animal shelters) and immune status. Transmission (direct contact, by contaminated food, water, objects, grooming and vectors like flies etc.) happens fecal-orally by ingestion of highly infectious and very resistant cysts being discharged by other animals or humans. Only five to ten cysts are enough to cause an infection.

G. duodenalis has an asexual life cycle. In the duodenum of the infected animals, two growth forms, so-called trophozoites, emerge from the incorporated cysts (excystment). These multiply by duplication and attach via suckers to the duodenal surface. Free trophozoites turn into their lasting forms, the cysts (encystment), especially in the ileum. These are excreted in large amounts (10⁷/g feces) and mostly intermittent, i.e. not with every defecation. The prepatent period averages 5 to 16 days.

The main symptom of giardiasis is diarrhoea, more or less intensive, that can run from symptomatic (acute, chronic, self-limiting, periodic-intermittent or continuous) to asymptomatic. Independent on the progression, cysts and/or trophozoites can be excreted (primarily with strong diarrhoea). *Giardia* cysts can be differentiated from cysts of different coccidia species only by microscopical experienced people. This is similarly true for *Giardia*- and *Tritrichomonas foetus* trophozoites.

Numerous studies show the advantages of the enzyme immunoassay for the qualitative detection of *Giardia duodenalis* compared with conventional detection methods for *Giardia* like the microscopic proof of cysts and/or trophozoites in the direct fecal smear or after enrichment by flotation or in the modified sedimentation ("MIFC method") or even via intestine biopsy.

Based on its high sensitivity and specificity, the **MegaELISA® GIARDIA** perfectly fits as standard test for veterinary laboratories for the fast and standardised routine test of intact *G. duodenalis* cysts, trophozoites and their disintegration products.

2. INTENDED USE

MegaELISA® GIARDIA is an enzyme immunoassay for the qualitative detection of *Giardia duodenalis* surface antigens in the feces of pocket pets, pets and farm animals.

3. TEST PRINCIPLE

In the **MegaELISA® GIARDIA**, specific antibodies are applied in a sandwich principle. Monoclonal antibodies against specific *Giardia* antigens are bound to the surface of the well in the ELISA plate. A suspension of the feces sample to be tested as well as internal controls are incubated in a first step at room temperature (20–25°C) with biotinylated anti-*Giardia* antibodies (Conjugate 1) in the wells of the ELISA plate. During presence of *Giardia* antigens, a sandwich complex results, consisting of immobilised antibodies, the *Giardia* antigen and the biotinylated antibody. After a wash step, Streptavidin poly-peroxidase conjugate (Conjugate 2) is added and incubated again at room temperature (20–25°C). Non-bound Streptavidin poly-peroxidase conjugate is removed in another wash step. After addition of Substrate, the bound enzyme changes the colourless solution in the wells of the ELISA plate with positive samples into a blue solution. By addition of Stop solution, the colour changes from blue to yellow. The extinction is proportional to the concentration of *Giardia* antigens present in the sample.

4. TEST-KIT COMPONENTS

4.1. REAGENTS / COMPONENTS PROVIDED WITH THE TEST-KIT

The reagents of one package are sufficient for 96 determinations.

	96's kit	Content
1.	12 separable strips	ELISA plate with separable microtitre strips in holding frame with 8 determinations each, coated with monoclonal antibodies against <i>G. duodenalis</i>
2.	100 ml	Sample dilution buffer (SDB, protein buffered NaCl solution, blue stained), ready to use, white cap
3.	100 ml	Wash buffer (10× phosphate buffered NaCl solution, contains 0.1 % thimerosal), brown cap
4.	2 ml	Positive Control (inactivated <i>G. duodenalis</i> antigen), ready to use, red cap
5.	2 ml	Negative Control , ready to use, blue stained, white cap
6.	13 ml	Conjugate 1 (monoclonal Biotin conjugated anti- <i>G. duodenalis</i> antibody in stabilised protein solution, red stained), ready to use, orange cap
7.	13 ml	Conjugate 2 (Streptavidin poly-peroxidase conjugate in stabilised protein solution), ready to use, orange stained, white cap
8.	13 ml	Substrate (3,3',5,5'-Tetramethylbenzidine [TMB]/H ₂ O ₂), ready to use, blue cap. Do not expose to sun light.
9.	12 ml	Stop solution (1 M Sulfuric acid), ready to use, yellow cap

	96's kit	Content
10.	1	Cover foil
11.	1	Evaluation sheet
12.	1	Instructions for use

4.2. REQUIRED, NON-PROVIDED REAGENTS / DEVICES

- Distilled or deionised water
- Disposable sample tubes
- Disposable plastic pipettes
- Laboratory pipette (range between 5 and 1000 μ l)
- Vortex mixer
- Filter paper/laboratory cloths
- Measuring cylinder (1000 ml)
- Stop watch
- Washing station for ELISA plates or multichannel pipette (300 μ l)
- Photometer for ELISA plates (450 nm, possibly reference filter 620)
- Waste bin with a 0.5% hypochlorite solution

5. STABILITY AND STORAGE

2-8°C

Store at
2-8°C

Expiry date
– see label

For veterinary use
only

Lot number

In vitro diagnosticum

Do not use test-kit
components from
different kits, lot
numbers or beyond
stated expiry date.

Follow instructions for
use precisely

- Microbial contamination must be avoided. After the expiration date, the guarantee of quality is not provided.
- When opening the sealed aluminium pouch with the ELISA plate, care must be taken that the reclosable clip closure is not damaged or even separated.
- A direct influence of light onto the colourless Substrate must be absolutely avoided to prevent degradation or colour change to blue by autooxidation. If the Substrate has changed to blue, it cannot be used any more (functional incapacity).

6. INFORMATION ON THE SAMPLE MATERIAL

- Fecal samples of pocket pets, pets and farm animals can be used.
- Keep in mind that the sample material, as well as all used test-kit components, should have reached room temperature at the time of application.
- For the test, the required amount of feces as described in issue 7.3/Sample preparation, is needed. The amount depends on the consistency of the sample.
- Care must be taken that the fecal samples are collected in adequate transport tubes. Contact of the feces sample with transport media, preservatives, animal sera, metal ions, oxidising agents or detergents should be avoided, because these can lead to interferences with the **MegaELISA® GIARDIA**. When rectal smears should be used, take care that enough feces sample material (minimum 100 mg) is present for the test procedure.
- Due to the normally inhomogeneous or nest-like dissemination of antigens in the feces, the specimen material has to be mixed up homogeneously (spatula, vortex-mixer) before sampling.
- Non-cooled (20–25°C), the sample should be tested within 4 hours. At 2–8°C, the sample can be stored up to 3 days, permanently at minimum –20°C. Thawing and freezing for several times must be avoided. Thawed samples are to be mixed homogeneously before dilution (optimally with vortex mixer).
- During investigations of larger animal groups (same stable/household), feces samples of clinical inapparent animals/contact animals also should be tested to identify asymptomatic carriers.

7. TEST PREPARATION

! Read instructions for use carefully before starting the test. For the reliability of the results, it is important to follow the instructions for use exactly as described. Carry out the test in the given order and without delay. Use clean disposable tips for each pipetting step.

7.1. GENERAL

- Remove the microtitre strips from the aluminium pouches after they have reached room temperature. To avoid loss by evaporation, it is advisable to cover or mask them. Immediately return non-used microtitre strips into the pouch, close the pouch and store it at 2–8°C.
- Mix the reagents well immediately before their use.
- After use, reagents should be stored at 2–8°C again.
- Microtitre strips that were once used cannot be used again.
- Damaged test-kit components must not be used.
- Direct contact of samples with the test-kit components has to be avoided (cross-contamination).
- Direct solar radiation during testing must be avoided.

7.2. PREPARATION OF WASH BUFFER 1:10

Mix 1 part of the Wash buffer concentrate with 9 parts of distilled water. When crystals are visible in the Wash buffer concentrate, dissolve these by warming the Wash buffer concentrate in a water bath at 37°C. The diluted buffer will be stable for 4 weeks when stored at 2–8°C.

7.3. SAMPLE PREPARATION / SAMPLE DILUTION 1:11

- Add 1 ml (1000 µl) **MegaELISA® GIARDIA** Sample dilution buffer (SDB) into an appropriate sample tube.
- Aspirate ca. 100 µl of liquid feces with a disposable plastic pipette or take ca. 50–100 mg of pulpy-compact feces with a spatula or disposable inoculation loop and add the sample into the SDB.
- Homogenise the sample-buffer mixture (SBM) well (repeated mixture with the pipette or vortex mixer). Leave the sample to rest for a short time to allow coarse feces particles to settle/sediment.
- The supernatant won hereof can be directly used for the test.
- **ATTENTION:** During the use of an ELISA washing station, the supernatant must be free of particles. To guarantee this, centrifuge the SBM at 2500 g for 5 minutes.
- Controls are ready to use and need not be diluted.

8. TEST PROCEDURE

- Mix sample material well before application.
- **It is recommended to test controls and samples in duplicates.**
- Before beginning, the position of the patient samples and the Controls in the microtitre strips must be determined on the Evaluation sheet supplied with the kit.

8.1. SAMPLING AND FIRST INCUBATION

- Place the required amount of test strips into the holding frame.
- Add **100 µl** each of
 - **Positive Control**
 - **Negative Control**
 - **diluted samples (SBM)**in the desired order into the particular cavity of the microtitre strip.
- Add **100 µl of Conjugate 1** (Biotin conjugated antibody) into each cavity and **mix** the suspensions by tipping softly at the plate rim.
- Cover the ELISA plate and incubate for **30 minutes at room temperature (20–25°C)**.

8.2. WASHING – CAREFULLY and STRICTLY FOLLOWING THE INSTRUCTIONS!

MANUAL WASHING

- Remove the incubation solution by emptying into a waste bin with hypochlorite. Disposal must be made according to official regulations.
- Beat the ELISA plate carefully several times onto filter paper/laboratory cloth.

- Wash 5× with 300 µl wash buffer per cavity. Shake plate with wash buffer slightly, then beat it.
- **Drain completely between each wash step, then beat carefully several times on a dry unused spot.**
- Before the last beating, take care that the reagent for the next pipetting step is prepared.

AUTOMATICAL WASHING

- **ATTENTION:** SBM still containing any particles should be manually removed from the cavities by shaking off (danger of blockage of the washing needles).
- Observe correct setting of the washing station (5×300 µl).
- To reach optimal washing results use at least 300 µl Wash buffer per cavity and wash step.
- Aspire liquid completely after each wash step.
- After the final wash step, beat the plate carefully but thoroughly on filter paper/laboratory cloth. There should not be any residual moisture in the cavities.

8.3. SECOND INCUBATION

- Add **100 µl Conjugate 2** (Streptavidin poly-peroxidase conjugate) into each cavity of the ELISA plate.
- Cover the ELISA plate and incubate for **15 minutes at room temperature (20–25°C)**.

8.4. WASHING

- Wash as described in issue 8.2.

8.5. THIRD INCUBATION

- Add **100 µl Substrate (TMB)** into each cavity of the ELISA plate.
- Cover the ELISA plate and incubate for exactly **15 minutes at room temperature (20–25°C) and in the dark.**
- Add **50 µl Stop solution** into each cavity, in the same order and at the same rate as for the TMB substrate solution.

8.6. MEASUREMENT

- After careful mixing (tipping softly at the plate rim), the extinction can be measured at 450 nm. An additional measurement at the reference wavelength 620 nm is recommended.
- Adjust the ELISA Plate Reader to zero against air – without ELISA plate.
- If duplicates or multiple determinations were done, calculate the mean **extinction values.**
- **ATTENTION:** Strongly positive samples can cause blackish precipitates of the Substrate. Pre-dilution of the sample with physiological NaCl solution, for example 1+1, is recommended. Then multiply the extinction values by 2.

9. TEST EVALUATION

9.1. REAGENT STABILITY – INTERNAL QUALITY CONTROL

During each testing, Positive and Negative Controls should be performed to ensure the stability of reagents as well as the correct test procedure.

The test was correct, when following reference values were reached:

- **Negative Control** (450 nm) **OD < 0.2** or (450/620 nm) **OD < 0.160**
and

- **Positive Control** (450 nm or 450/620 nm) **OD > 0.8**

If the OD of the negative control is 450 nm > 0.2 (450/620 nm OD > 0.160), the washing process was probably insufficient.

A variation of the required values as well as any clouding or blue staining of the colourless Substrate before adding into the cavities can be a hint on expiry of reagents. If the prescribed issues are not fulfilled, following should be ensured before test repetition:

- Stability of the used reagents
- Functionality of the used machines (calibration)
- Visual control of the test-kit components for contamination or leakage; a blueish stained Substrate solution should not be used.

9.2. CUT-OFF

For determination of the Cut-off, 0.15 extinction units are added to the measured extinction of the Negative Control.

Cut-off = Extinction of Negative Control + 0.15

9.3. EVALUATION OF TEST RESULTS

Positive test result

Sample with an extinction more than 10 % above the calculated Cut-off.

Marginal test result

Sample with an extinction in the field of 10 % above and below the calculated Cut-off. A repetitive test with a new sample being marginal, too, should be evaluated as negative.

Negative test result

Sample with an extinction more than 10 % below the calculated Cut-off.

9.4. INTERPRETATION OF THE TEST RESULTS

- The diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptoms as well as serological data.
- Bacterial contamination or repeated freeze-thaw cycles of the specimen may affect the extinction values.
- The **MegaELISA® GIARDIA** detects antigens of *Giardia duodenalis* in the feces of pocket pets, pets and farm animals. A relationship between the magnitude of a detected value and the onset or severity of clinical symptoms can not be deduced therefrom.

- A positive result does not exclude the presence of other infectious agents.
- A negative result does not exclude a possible *Giardia* infection. In case of existing diarrhoea problem and one-time negative test, a new test with a new feces sample (optimally collection feces sample: individual testing of at least three consecutive feces samples) should be carried out within 7 days.
- A marginal test result can be caused by an inhomogeneous distribution of the *Giardia* in the feces sample. In this case, a second suspension from the same feces sample should be tested, or a further feces sample should be requested for examination.

10. SAFETY MEASURES AND PRECAUTIONS

- Wear safety clothes (disposable gloves, lab coat, safety glasses) and wash the hands after finishing the test procedure.
 - No eating, drinking and smoking during test procedure.
 - The sample material and Controls must be seen as potentially infectious and disposed of accordingly after test procedure, together with the used test-kit components.
 - Do not use reagents of other manufacturers together with the reagents of this test-kit.
 - Do not mix reagents of different Lot numbers.
 - Do not change the caps of the individual reagents among each other.
 - Close bottles immediately after use to avoid evaporation and microbial contamination.
 - After first opening, the reagent bottles must be checked for microbial contamination before reusing.
 - Wash buffer contains Thimerosal as preservative.
 - Stop solution contains 1 M Sulfuric acid. Sulfuric acid irritates eyes and skin. Upon contact with the eyes, rinse thoroughly with water and consult a doctor. Keep out of the reach of children.
 - Substrate solution contains Hydrogen peroxide.
- Absolutely avoid eye, skin or mucosa contact.**
- **ATTENTION:** Liquid waste containing Stop solution must be neutralised prior to addition to a hypochlorite solution.
 - Only use clean pipette tips, dispenser and laboratory material.
 - Never pipet with your mouth.
 - Label sample material and associated sample tube to ensure a precise assignment.
 - Use a new sample tube for each sample.
 - During the incubation and wash steps, the cavities must be completely covered with liquid.
 - The cavities of the ELISA plate must not stay empty longer as necessarily required.

- To avoid cross contamination and falsely elevated rest result, pipet the patient samples and Conjugate carefully into the cavities of the ELISA plate.
- **MegaELISA® GIARDIA** is only provided for the use of qualified staff that knows the work technique objectionable.

The test procedure, the information, the safety measures and warnings in the Instructions for use are to be followed strictly. For use of the test-kit in diagnostical instruments, the test method has to be validated. Each change in appearance, composition and the test procedure as well as each use in combination with other products not authorised by the manufacturer is not permitted; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents being associated herewith. No liability is accepted for false test results due to visual evaluation.

11. LIABILITY

The entire risk of liability in connection with the use of this product remains with the purchaser. The manufacturer is not liable for indirect, special or consequential damages of any kind resulting from the usage, test performance and evaluation of this product.

12. SHORT PROTOCOL MegaELISA® GIARDIA

TEST PREPARATION

Warm needed test strips and reagents to **room temperature (20–25 °C)**.

Dilute Wash buffer 1:10 with distilled water.

Dilute feces samples 1:11 with Sample dilution buffer (100 µl feces + 1000 µl SDB)

TEST PROCEDURE

Add **100 µl of Positive and Negative Controls**
or **1:11 diluted feces samples**
each in a cavity.

▼
Add **100 µl of Conjugate 1**
in each cavity. Mix carefully.

▼
Cover and incubate **30 minutes at room temperature (20–25 °C)**.

▼
Wash 5× with soft shaking with **300 µl Wash buffer each**.

▼
Add **100 µl of Conjugate 2**
in each cavity. Mix carefully.

▼
Cover and incubate **15 minutes at room temperature (20–25 °C)**.

▼
Wash 5× with soft shaking with **300 µl Wash buffer each**.

▼
Add **100 µl of Substrate**
in each cavity. Mix carefully.

▼
Cover and incubate **15 minutes at room temperature (20–25 °C)**
in the dark.

▼
Add **50 µl of Stop solution**
in each cavity.

▼
Photometric evaluation of the cavities at 450 nm.

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