




# MegaELISA® VET

**Veterinary tests with laboratory precision**

**Enzyme immunoassay for the qualitative detection of specific IgG antibodies, antigens or proteins**



- **Standardised test procedure**
- **Manual or automated processing / evaluation**
- **Parallel analysis of large sample numbers**  
with plasma or serum / with feces
- **Routine diagnostics**
- **Format: 48 or 96 wells; ready to use**

## Test principle


Antibody detection  
(*Leishmania*, *Anaplasma*)




Antigen detection  
(*Giardia*, *Cryptosporidium*)



## Test procedure



## Test evaluation



Art. No.

600048EK1	<b>MegaELISA® ANAPLASMA (48 tests)</b>
600096EK1	<b>MegaELISA® ANAPLASMA (96 tests)</b> Qualitative detection of IgG antibodies against <i>Anaplasma phagocytophilum</i> in plasma or serum of the dog
870096EG1	<b>MegaELISA® CRYPTO (96 tests)</b> Qualitative detection of <i>Cryptosporidium parvum</i> in feces of pocket pets, pets and farm animals
740096EG1	<b>MegaELISA® CRYPTO-GIARDIA (96 tests)</b> Qualitative detection of <i>Cryptosporidium parvum</i> and <i>Giardia duodenalis</i> in feces of pocket pets, pets and farm animals
991048EK1	<b>MegaELISA® EHRLICHIA canis (48 tests)</b>
991096EK1	<b>MegaELISA® EHRLICHIA canis (96 tests)</b> Qualitative detection of IgG antibodies against <i>Ehrlichia canis</i> in plasma or serum of the dog
821048EK1	<b>MegaELISA® FCoV (48 tests)</b>
821096EK1	<b>MegaELISA® FCoV (96 tests)</b> Qualitative detection of IgG antibodies against the Feline Coronavirus (FCoV) in plasma, serum or effusion of the cat
850096EG1	<b>MegaELISA® GIARDIA (96 tests)</b> Qualitative detection of <i>Giardia duodenalis</i> in feces of pocket pets, pets and farm animals
987048EK1	<b>MegaELISA® LEISH canine (48 tests)</b>
987096EK1	<b>MegaELISA® LEISH canine (96 tests)</b> Qualitative detection of antibodies against <i>Leishmania infantum</i> in plasma or serum of the dog

### **Distribution products**

---

**Bionote**

Anigen AIV Ab, Anigen PRRS Ab

**TECO**

CRP, Haptoglobin

**MEGACOR Diagnostik GmbH**

Lochauer Str. 2 – 6912 Hörbranz – AUSTRIA

[www.megacor.com](http://www.megacor.com) – [info@megacor.at](mailto:info@megacor.at) – [+43 5573 85400](tel:+43557385400) – [+43 5573 85400-4](tel:+435573854004)