




# MegaELISA® VET

**Veterinary tests with laboratory precision**


**Enzyme immunoassay for the qualitative detection of specific IgG antibodies, antigens or proteins**




- **Standardised test procedure**
- **Manual or automated processing / evaluation**
- **Parallel analysis of large sample numbers**  
with plasma or serum / with feces
- **Routine diagnostics**
- **Format: 48 or 96 wells; ready to use**

## Test principle


Antibody detection  
(*Leishmania*, *Anaplasma*)




Antigen detection  
(*Giardia*, *Cryptosporidium*)



## Test procedure



## Test evaluation



Art. No.

600048EK1	<b>MegaELISA® ANAPLASMA</b> phagocytophilum (48 tests)
600096EK1	<b>MegaELISA® ANAPLASMA</b> phagocytophilum (96 tests) Qualitative detection of IgG antibodies against <i>Anaplasma phagocytophilum</i> in plasma or serum of the dog
760096EE1	<b>MegaELISA® CRP</b> canine (96 tests) Quantitative detection of C-reactive Protein (CRP) in plasma or serum of the dog
870096EG1	<b>MegaELISA® CRYPTO</b> (96 tests) Qualitative detection of <i>Cryptosporidium parvum</i> in feces of pocket pets, pets and farm animals
989096EG1	<b>MegaELISA® DIRO</b> Antigen (96 tests) Qualitative detection of <i>Dirofilaria immitis</i> antigens in plasma or serum of the dog
991048EK1	<b>MegaELISA® EHRLICHIA</b> canis (48 tests)
991096EK1	<b>MegaELISA® EHRLICHIA</b> canis (96 tests) Qualitative detection of IgG antibodies against <i>Ehrlichia canis</i> in plasma or serum of the dog
825048EK1	<b>MegaELISA® ENCEPHALITOOZON</b> cuniculi (48 tests)
825096EK1	<b>MegaELISA® ENCEPHALITOOZON</b> cuniculi (96 tests) Qualitative detection of IgG antibodies against <i>Encephalitozoon cuniculi</i> in plasma or serum of the rabbit
821048EK1	<b>MegaELISA® FCoV</b> (48 tests)
821096EK1	<b>MegaELISA® FCoV</b> (96 tests) Qualitative detection of IgG antibodies against the Feline Coronavirus (FCoV) in plasma, serum or effusion of the cat
950048EG1	<b>MegaELISA® FeLV</b> Antigen (48 tests)
950096EG1	<b>MegaELISA® FeLV</b> Antigen (96 tests) Qualitative detection of Feline Leukaemia Virus (FeLV) antigens in plasma or serum of the cat
804096EK1	<b>MegaELISA® FIV</b> (96 Tests) Qualitative detection of antibodies against the Feline Immunodeficiency Virus (FIV) in plasma or serum of the cat
850096EG1	<b>MegaELISA® GIARDIA</b> (96 tests) Qualitative detection of <i>Giardia duodenalis</i> in feces of pocket pets, pets and farm animals
987048EK1	<b>MegaELISA® LEISH</b> canine (48 tests)
987096EK1	<b>MegaELISA® LEISH</b> canine (96 tests) Qualitative detection of antibodies against <i>Leishmania infantum</i> in plasma or serum of the dog

## **Distribution products**

---

### **Bionote**

AIV Ab ELISA  
PRRS Ab ELISA

### **TECO**

Canine CRP ELISA  
Haptoglobin Equine ELISA  
Haptoglobin Feline/Canine ELISA

## **MEGACOR Diagnostik GmbH**

Lochauer Str. 2 – 6912 Hörbranz – AUSTRIA

[www.megacor.com](http://www.megacor.com) – [✉ elisa@megacor.at](mailto:elisa@megacor.at) – [☎ +43 5573 85400](tel:+43557385400) – [📠 +43 5573 85400-4](tel:+435573854004)