

POCKIT[®] VET ad us. vet. – Who's who

Sample preparation: DNA/RNA extraction

1. Automatic extraction with magnetic beads

taco[™]



taco[™]

- 24 samples at a time
- 372001RG1



taco[™] mini

- 8 samples at a time
- 381001RG1

PetNAD[®]

2. Manual extraction

- PetNAD[®] extraction kit 375050PX2
- Grind-N-Go 348050PS2
- IQ Plus[™] 360050PX1



POCKIT[™] Central

- extraction and amplification in one device
 - total duration about 90 minutes
- 488001PG1

Test procedure: DNA amplification

POCKIT[™] PCR



POCKIT[™] Micro Duo

- up to 4 samples
 - evaluation using two wavelengths
 - Micro SD card
 - duration about 45 minutes
- 389001PG1



POCKIT[™] Micro Plus

- up to 4 samples
 - duration about 45 minutes
- 383001PG1

- ready-to-use reagents



For detailed list of kits for pathogen identification please visit GeneReach website

→ https://www.genereach.com/index.php?func=product&action=view&product_no=9 (Scroll to "Download")

Test evaluation



qualitative: Yes/No